



Parameter	Model	INS-FH800EP	INS-FH1200	INS-FH1500	INS-FH1800
Wind speed		0.25-0.45m's range adjustable			
Intake pressure		>0.5Pa			
Power supply		AC single phase 220V/50HZ			
Maximum power		400W	400W	800W	1000W
Weight		100KG	130KG	180KG	1330KG
Working size (mm)W1×D1×H1		600*500*500	988*655*1200	1288*656*1200	1588*665*1200
Overall size (mm)/W-D-H		760-580*1770	1200*850*2350	1500*850*2350	1800-850*2350
UV lamp specification and quantity		28W*1	28W 1	30W 1	40W 1