

Clean Laminar Flow Booth



Parameter	Model	INS-CL-700	INS-CL-1300	INS-CL-1900	INS-CL-2500
Purification area (mm)		650*610	1250*950	1850*1220	2450*1220
Overall dimension (mm)		700*750*750	1300*1100*750	1900*1400*750	2500*1400*750
Power consumption		0.4	0.8	1.2	1.6
Gel-type HEPA filter		700*500*90mm*1	1000*500*90mm*2	1000*500*90mm*3	1000*550*90mm*4
Illuminance		300LUX	300LUX	300LUX	300LUX
DOP test port		YES	YES	YES	YES
Diferential pressure gauge		YES	YES	YES	YES