

Clean Sampling Vehicle



Parameter	Model	INS-JC-900	INS-JC-1200
Cleanliness		Level 100 @ $\geq 0.5 \mu\text{M}$ (Federal 209E)	
Number of colony		≤ 0.5 pieces/dish hour (90mm culture plate)	
Average Wind speed		0.25 0.45m/s	
Noise		$\leq 62\text{dB}$ (A)	
Vibration half peak		$\leq 5\mu\text{M}$ (x, y, z)	
Illuminance		$\geq 300\text{LX}$	
Power supply		AC single phase 220V/50HZ	
Maximum power consumption		800W	
Weight		< 3	
Internal size (mm)W1×D1×H1		730*635*930 (1300)	1030*635*930 (1300)
External size (mm)WxDxH		885*690*1910	1185*690*1910
HEPA filter specifications and quantity		725*600*50*1	1025*600*50*1
LED/UV lamp specifications and quantity		20W*1/20W*1	30W*1/30W*1