

Co2 Incubator



Model	INS-CO-140	INS-CO-140	INS-CO-140	INS-CO-140
Type	Air jacket		Water jacket	
Capacity	80L	160L	80L	160L
Temp. range	RT+5~50°C			
Temp. fluctuation	±0.2°C			
Temp. resolution	0.1°C			
CO2 range	0~20%			
Humidification mode	Natural evaporation			
Voltage	AC 220V 50HZ			
Input power	450w	720W	680W	950W
Internal size mm	400×400×500	500×500×650	400×400×500	500×500×650
External size mm	530×505×785	630×605×935	530×505×785	630×605×935