

# HandScan Ultrasonic Scanner



<b>Model</b>	<b>INS-0185</b>
Weight	701 g (including battery)
Device Type	Handheld ultrasound scanner
Power	Battery and AC powered, >3 hours continuous
Probe	Multi-frequency waterproof mechanical sector probe
Protection	Rubber protector against dirt and impact
Operation	Shoulder and waist belts for hands-free use
Durability	Waterproof and dustproof
Standard Package	Main unit, battery, AC adapter, probe, rubber protector, shoulder/waist belts, carrying case
Optional Accessories	Multi-language support