



Parameter \ Model	INS-CJ-1G	INS-CJ-2G	INS-CJ-1B	INS-CJ-1C	INS-CJ-1BU (OPEN TYPE)	INS-CJ-1CU (OPEN TYPE)
Cleanliness level	Level 100 @≥ 0.5 μM (Federal 209E)					
Colony count	≤0.5 pieces/dish hour (90mm culture plate)					
Average wind speed	0.25 0.45m/s					
Noise	≤62dB (A)					
Vibration half peak	≤5μM (x, y, z direction)					
Illuminance	≥300LX					
Power supply	AC single phase 220V/50Hz					
Cleanliness level	100W	500W	490W	640W	490W	640W
Colony count	80KG	110KG	90KG	140KG	90KG	160KG
Inner size (WxDxH) mm	700x410x580	1150x410x580	880x510x580	1370x510x580	900x510x600	1770x510x600
External size (WxDxH) mm	860x590x1600	1310x590x1600	1040x 690x1600	1530x690x1600	930x740x1440	1800x510x1440
Applicable number of people	Single Person/ Single Side	Double Person/ Single Side	Single Person/ Single Side	Double Person/ Single Side	Single Person/ Single Side	Double Person/ Single Side
Specification and quantity of high-efficiency filters	680x565x50	1130x565x50	880x510x50	1350x565x50	830x560x50	1700x560x50
Specification and quantity of fluorescent/ultraviolet lamps	12W x i20W x 1	12W x i30W x 1	12W x i20W x 1	20W x i36W x 1	12W x i20W x 1	20W x i36W x 1