

IL-35-05 Technical Table

Parameter	Specification
Power Supply	380V, 50Hz
Rated Power (kW)	1.0
Temperature Range (°C)	50 – 300
Temperature Fluctuation (±°C)	±1
Temperature Deviation (°C)	+5
Work Room Size (mm W×D×H)	350 × 350 × 350
Overall Dimensions (mm W×D×H)	800 × 800 × 1000
Inner Material	Stainless steel
Notes	Microcomputer temperature control, hot air circulation, toughened glass viewing window

