

IL-70-05 Technical Table

Parameter	Specification
Power Supply	380V, 50Hz
Rated Power (kW)	2.0
Temperature Range (°C)	50 – 300
Temperature Fluctuation (±°C)	±1
Temperature Deviation (°C)	+5
Work Room Size (mm W×D×H)	450 × 550 × 550
Overall Dimensions (mm W×D×H)	1000 × 850 × 1200
Inner Material	Stainless steel
Notes	Microcomputer temperature control, hot air circulation, toughened glass viewing window

