

INS-CL-1900 Technical Table

Parameter	Specification
Cleanliness level	Class 100 (US Federal Standard 209E)
Average wind speed	0.25 – 0.45 m/s
Noise level	≤ 62 dB(A)
Power supply	AC 220V / 50Hz
Purification area (mm)	1850 × 1220
Overall dimensions (mm)	1900 × 1400 × 750
Power consumption	~1.2 kW
HEPA filter	Gel-type HEPA, 1000 × 500 × 90 mm × 3
Illuminance	≥ 300 Lux
DOP test port	Yes
Differential pressure gauge	Yes

