

Illumination Incubator



Model	LC-70-I	LC-150-I	LC-250-I	LC-360-I	LC-460-I
Capacity	70L	150L	250L	360L	460L
Constant temperature method	Forced convection, heating and cooling balance				
Temp. range	10-50 when illuminated 0-50 °C when not illuminated				
Temp. resolution	0.1°C				
Temp. fluctuation	0.5°C				
Temp. accuracy	±1°C (at 37°C)				
Light source illuminance	0-12000LX (LED light source)				
Internal material	Stainless steel				
External material	cold rolled steel plate is coated with antibiotic resistant coating				
Heater	Stainless steel heater				
Display	LCD display				
Voltage	AC 220V 50HZ				
Power consumption	0.4KW	0.6KW	0.8KW	0.9KW	1.0KW
Internal size (DxWxH)mm	350×400×500	400×500×750	500×550×900	600×600×1000	650×650×1100
External size (DxWxH)mm	540×540×1070	640×640×1360	720×690×1580	820×740×1680	870×790×1780