

Integrated Muffle Furnace



Model	Related Power (Kw)	Voltage	Related Temperature (C)	Size of Workroom (mm)
INS-2.5-10	2.5	220	1000	200*120*80
INS-4-10	4	220	1000	300*200*120
INS-8-10	8	380	1000	400*250*160
INS-12-10	12	380	1000	500*300*200
INS-2.5-12	2.5	220	1200	200*120*80
INS-5-12	5	220	1200	300*200*1200
INS-10-12	10	380	1200	400*250*160
INS-12-12	12	380	1200	500*300*200
INS-4-13	4	220	1300	250*150*200
INS-6-13	6	380	1300	250*150*100
INS-8-13	8	380	1300	500*200*180
INS-10-13	10	380	1300	400*200*160