

Vertical Autoclave



Model	Volume	Power Supply	Power	Dimension of sterilizing chamber	Weight	Package size
INS-IL-30-05	30L	220V	3.5KW	φ280*400MM	47Kg	510*510*1170MM
INS-IL-40-05	40L	220V	3.5KW	φ340*425MM	55Kg	570*570*1250MM
INS-IL-50-05	50L	220V	3.5KW	φ340*525MM	58Kg	570*570*1300MM
INS-IL-70-05	75L	380V	4.5KW	φ380*530MM	77Kg	630*630*1400MM
INS-IL-100-05	100L	380V	6.5KW	φ500*560MM	120KG	660*630*1420MM
INS-IL-150-05	150L	380V	6.5KW	φ500*760MM	130KG	660*66*1200MM
INS-IL-200-05	200L	380V	8.5KW	φ600*760MM	180Kg	790*790*1480MM