

Vertical Blast Drying Oven



Model (Stainless Steel Inner)	Power	Temperature Range	Temperature fluctuation	Temperature fluctuation	Consumed Power	Dimensions (mm)
INS-IL-35-05	220V 50HZ	50°C~300°C	±1°C	±5°C	1000W	350*350*350
INS-IL-50-05					1500W	350*450*450
INS-IL-70-05					2000W	450*550*550
INS-IL-100-05	2500W				500*600*750	
INS-IL-150-05	380V 50HZ				4500W	800*800*1000