

# Distillation Unit

Cat. No : LS-21007



## FEATURES

- Horizontal thick wall heat resistant borosilicate glass boiler fitted with **STAINLESS STEEL HEATER**.
- High efficiency coil condenser ensure low temperature, high output quality distillate of conductivity 3 to 4µs/cm.
- Compact construction suitable for bench and wall mounting.
- Screwthread connectors enable easy fitting or removal of hose safely.
- Built-in constant level device with Vendura-Nova valve for drainage provides easy descaling without dismantling.
- Safety thermostat built into the heater protect the Favorit water still in the event of water supply failure and automatically resets when the water supply is restored.

## SPECIFICATIONS

<b>Output: litres/hour</b>	4L/hour	<b>Water Supply</b>	1 litre / min Min. pressure 3 psi (0.2 x 10 <sup>5</sup> NM <sup>-2</sup> )
<b>Quality of distillate:</b>		<b>Power</b>	220 / 240V 50 / 60 Hz single phase 13A fused
pH	5.0 - 6.5	<b>Dimensions :</b> (w x d x h)	50 x 15 x 45cm
Conductivity µS/cm	3.0 - 4.0	<b>Catalogue Number:</b>	LS-21007
Resistivity megOhm-cm	0.25 - 0.30		
Pyrogen content	Pyrogen free		
Temperature	35°C-50°C		